
Fluorescence Observation Kit

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
Nomenclature / 術語表

- ① UV/ Blue light Shield is designed with trigger function for protection to avoid the leak of blue light to damage eyes.
遮光擋板具有起動保護功能設計，以避免洩漏的藍光，對眼睛造成損害。
【 Note 】 Starting and operation requires keeping insert the light shield; if draw it out, lighting will be off immediately;
Reset need to turn on lighting again.
【 註 】 啟動和運行，需要保持插入遮光擋板，若取出，照明將立即關閉，重新插入後需再次開啟照明。
- ② Emission (Barrier) filter block shorter wavelengths and have high transmission for longer wavelengths.
The extention tubes (can be freely increased or decreased, up to 4) enable the filter close to the object lens then reduce glare.
(Two kind of length sized tube prepared for different size of microscope stand arms)
發射(屏障)濾鏡阻斷較短的波長，僅傳輸較長的波長。增距接環(可自由增減，至多可使用到4個)，使濾鏡靠近物鏡，減少眩光。
- ③ Slide control: 5-steps lighting intensity, and memory function for last time setting.
滑動控制：5段照明強度，最後一次設置記憶功能。
- ④ Illumination angle is adjustable. Turn both knobs and adjust different angles, to find out the proper lighting angle for observation.
照明角度可調整。同時轉動旋鈕，調整不同的角度觀察，以找出適當的照明角度。
【 Note 】 When the lighting angle is 45 degree, the lighting effect is best but glare are also strong; If turn it to 0 degree, then the lighting is almost paralleled to chamber slide, which usually helps to avoid the reflective glare.
【 註 】 當照明角度45度，照明效果是最好的，但眩光也強，如果為0，那麼光線幾乎平行於觀察玻片皿，這通常有助於避免反射眩光。
- ⑤ Chamber Slide: 4 kinds of chamber slides are available for different purposes.
(2 types standards : 3 large/medium/small wells and 10 small wells)
觀察玻片皿：共有四種型式供選擇使用(其中大/中/小三孔和十小孔為標準配備)。
- ⑥ Manual iris door 手動光量控制門
【 注意 】 光量控制以不遮住物鏡視野，造成暗角的情況下請。

Precautions for use 使用注意事項

⚠ Warning 警告 ⚠

Please ensure that you read and understand the safety precautions before using.

To avoid the fire, electric shocks or result in a body injury, please obey the following instructions.

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1. Keep this device from the reach of children.

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2. This device equipped high-brightness excitation LEDs; it is intense enough to injure eyes. While operating, for eyes safety, DO NOT look directly excitation light source, or through the slot of device to look into, or stare the leak of light at close range.

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3. DO NOT leave this device near the heat source or use in a hot environment; DO NOT block the ventilation hole / vent hole of cooling fan. This device has ventilation holes at the front/ rear / bottom to prevent the internal temperature from becoming too high. If the inner temperature rises too high (over 70°C / 158°F), the lighting will be shut down for the over-heating protection until the temperature cools down and falls below the set point temperature (about 1 minute), it will automatically open the circuit and power supply. If over-heating happened repeatedly, please contact your distributor.

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4. Although it's allowed to place the water drops in chamber slide for using, and the device has power circuit breaker designed to prevent large amounts of droplets may cause device damage or danger; you should still keep the device away from water or other liquids. A humid device could result in fire or electric shock. If it's wet, immediately turn off the device power, or unplug from outlet, then wipe out the liquid the device surface.

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
5. Turn off the LED lights when you aren't using them. 不使用LED燈時請關閉。

6. Pay attention to the running of cooling fan. 要注意冷卻風扇的運行。

7. Unplug this device when unused for long periods of time. 長時間不使用該裝置時，請拔掉插頭。

8. Refer all servicing to your distributor. Servicing is required when the device has been damaged in any way, such as power-supply cord or plug is damaged, liquid has been spilled, the apparatus has been exposed to rain or moisture, does not operate normally, or has been dropped. Do NOT disassemble, alter or modify to this device.

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Disclaimer

- Thank you for using this product. While every effort has been made to ensure that the information contained in this guide is accurate and complete, no liability can be accepted for any errors or omissions.
- We (Lumos Tech.) reserves the right to change the specifications described herein at anytime without prior notice.
- No part of this guide may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any language in any form, by any means, without the prior written permission of Lumos Technology Co., Ltd.

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In case of any question, please contact your
X-Loupe® distributor, or e-mail : service@x-loupe.com

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